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ENGLISH

IN ORDER TO ENSURE THAT YOU CAN USE THEM CORRECTLY PLEASE CAREFULLY READ THIS INSTRUCTION MANUAL BEFORE USING THE SCALES AND FOLLOW THE INSTRUCTIONS DURING USE.

Specifications

Capacity: 180 kg

BMI

Water
(10 to 75%, 0,1%)

Muscle
(0,1%)

Bonne
(0,1%)

Fat
(4 to 80%, 0,1%)

- *Equipped with a high precise strain gauge sensor system (4load cells)
- * Bamboo platform
- * LCD display size: 80 x 45 mm
- * Division: 100g/0.2lb
- * Tap to switch on
- * Individual data memory:10 people
- * Height range: 100-240cm
- * Age range: 10-80 years
- * Automatic zero resetting
- * Automatic switch off
- * Overload indicator
- * Low battery indicator
- * Power: 1 x CR2032 / CR2034 - 3V lithium battery
- * Recommended working environment: Temperature: 0°C -40°C/32°F-104°F; Relative Humidity <85%
- * Choice of weight unit: kilogram (kg); stone(st); pound (lb)

THIS SCALE WAS DESIGNED FOR HOUSEHOLD USE ONLY. IT MUST NOT BE USED FOR COMMERCIAL PURPOSES.

1. Remove the insulation slip from the battery contacts before use.
2. Select the weighing mode using the switch (kg) located on the back of the scales.
3. Place the scales on a firm, flat and even surface to ensure accuracy.
4. Tap the scales' platform to switch on the product. When the LCD display reads 0.0, step onto the platform, stand motionlessly while the scales compute your weight.
5. Your weight will be displayed and will flash twice before it is fixed.
6. When you step off the scales they will automatically switch off after 6 seconds.
7. If no weighing is carried out while 0.0 is displayed the scales will automatically switch off after 6 seconds.

Body Fat Measurement

In order to ensure accurate results, users should enter the required personal parameters as below: age, sex and height. These scales can store different personal parameters for up to 10 people. The technology used in these scales, Bio-electrical Impedance Analysis (BIA), is more technologically advanced than the older Body Mass Index (BMI) method in which fat% was calculated by using height and weight only. Our new BIA technology

uses 5 parameters to get a more accurate resolution for the analysis: weight and biological resistance, which are automatically measured by the scale, together with one set of personal parameters (age, gender and height), which are entered manually.

Mode selection and parameter setting:

Press and hold "SET" to enter parameter settings:

Press ▲ or ▼ to switch between Body parameters 0-9. Press "SET" to confirm your setting.

Press ▲ or ▼ to select male or female. Press "SET" to confirm your setting.

Press ▲ or ▼ to set your age. Press "SET" to confirm your setting.

Press ▲ or ▼ to set your height. Press "SET" to confirm your setting.

The values are flashed and recorded after you confirm the parameter settings.

Stand on the scales after 0.0 appears on the LCD. First the weight is recorded, then they display fat and water, muscle, bone, calories and BMI, indicating that the scales are analysing. Please do not move while this is happening. Once finished, the weight and body fat ratio, body water ratio, body muscle ratio, body bone ratio, calorie and BMI data will be displayed alternately 2 times before the scales automatically turn off. Compare the body fat reading against the "Body Fat Ratio Chart" to confirm your body health status.

Body Fat Ratio Chart (%)

Age	FEMALE (%)				MALE (%)			
	Low	Healthy	High	Obese	Low	Healthy	High	Obese
6 - 39	< 21	21 - 33	33 - 39	> 39	< 8	8 - 19	19 - 25	> 25
40 - 59	< 23	23 - 35	35 - 40	> 40	< 11	11 - 22	22 - 28	> 28
60 - 100	< 24	24 - 36	36 - 42	> 42	< 13	13 - 25	25 - 30	> 30

Due to the natural differences between people these charts should be used for reference purposes only.

Water Ratio Chart & Muscle Ratio Chart (%)

Age	FEMALE (%)				MALE (%)			
	Water			Body Muscle	Water			Body Muscle
	Low	Healthy	High		Low	Healthy	High	
6 - 15	< 57	57 - 67	> 67	> 34	< 58	58 - 72	> 72	> 40
16 - 30	< 47	47 - 57	> 57		< 53	43 - 67	> 67	
31 - 60	< 42	42 - 52	> 52		< 47	47 - 61	> 61	
31 - 60	< 37	37 - 47	> 47		< 42	42 - 56	> 56	

Due to the natural differences between people these charts should be used for reference purposes only.

Bone Weight Chart

	FEMALE			MALE		
Weight	< 45Kg	45 - 60Kg	> 60Kg	< 60Kg	60 - 75Kg	> 75Kg
Bone weight	-5,8Kg	-6,2Kg	-6,5Kg	-6,5Kg	-6,9Kg	-7,2Kg

Due to the natural differences between people these charts should be used for reference purposes only.

Trouble shooting

1. When the LCD screen displays LO, you should install 1 (one) new CR2032 battery.
2. When the LCD screen displays Err or ERR1, it indicates the wrong procedure has been followed.
3. When the LCD screen displays OL, it indicates that the maximum weight is exceeded.
4. When the LCD screen displays ErrH, it indicates that fat measurement is over 80%.
5. When the LCD screen displays Err 2%, it indicates the wrong procedure has been followed for measuring body fat percentage.
6. When the LCD screen displays Err1, it indicates the fat measurement is below 4%.

Battery Safety Advice

ALWAYS:

- * Use the recommended batteries.
- * Make sure that the battery compartment is secure.
- * Fit batteries correctly, observing the plus and minus marks on the battery and compartment.
- * Store unused batteries in their packaging and keep away from metal objects that could cause a short circuit.
- * Remove dead batteries from the scales.
- * Remove batteries from the scales if they will not be used for long time.
- * Keep loose batteries away from young children.
- * Seek medical advice if a battery has been accidentally swallowed by young child.

NEVER:

- * Dispose of batteries in fires.
- * Attempt to recharge ordinary batteries.

Advice for use and care

- * Ensure the scales are put on a flat, even and firm surface. Avoid carpet or soft surfaces.
- * Always weigh yourself with the same scales placed on the same surface. Uneven floors may affect the reading.
- * The surface of the scales will be slippery if wet. So please ensure they are dry to ensure your safety when weighing!
- * Weigh yourself without shoes and socks, before meals and always at the same time of each day.
- * For better accuracy, it is not advised to take measurements within 2 hours of waking up.
- * Keep the scales in a dry place to protect their electronic components.
- * Clean after use with a soft damp cloth, do not use solvents or immerse the product in water.
- * Avoid overloading, otherwise, permanent damage may occur.
- * Do not disassemble the scales other than to replace the battery.
- * Avoid excessive impact or vibration to the unit, such as dropping the scales onto the floor, or dropping other objects onto the scales.
- * A user who has any metal material inserted into his/her body may get differences in body fat measurement readings.
- * People with pacemakers are advised not to use this device.

Changing the battery

1. Open the battery box lid on the rear face of the scales by releasing the clip and lift off the box lid.
2. Remove the used batteries and replace with a new set. Ensure that the polarities are correct.
3. After replacing the battery, close the battery compartment lid.

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